All SPOT RT CMOS cameras are compatible with:





Take a closer look.



New SPOT RT sCMOS Camera is a game changer.

The SPOT RT sCMOS uses Sony's breakthrough Pregius[™] CMOS sensor. Now you can experience unprecedented speed and sensitivity in a CMOS camera. Deep cooling allows dim images to be seen without becoming obscured by dark current. A global shutter ensures undistorted images of moving specimens. The RT sCMOS is optimal for fluorescence microscopy, FISH, GFP imaging, immunofluorescence and 3D deconvolution applications.

Cooling

Deep-cooling of the image sensor to -25°C supports long exposure times. *How cool is that*?

- Lower dark current
- Long exposure capability
- Software enable/disable cooler and fan

Sensitivity

Step up to a new level of sensitivity: 8x more sensitive than the previous RT3.

- Back-Illuminated sensor increases sensitivity
- 1/8 the read noise
- Enhanced Near-IR



Resolution

2 and 1/2 times the resolution of the RT3.

- 5 MP live and captured image resolution
- High resolution picture quality

Speed

Our RT3 camera was fast. This one's over seven times as fast. *Why wait?*

- Super fast USB3
- 36 frames per second at full resolution
- 70 frames per second with 2 x 2 binning

The suggested applications for this camera include:

Cell Biology

- Spinning Disk Confocal
- High Speed Multicolor Fluorescence
- FRAP

Intrinsic Imaging

- Ion Transport Physiology
- Electrophysiology
- Calcium Imaging
- Ratiometrics Imaging
- Voltage Sensitive Dyes



SPOT RT Cooled Camera with Sony Pregius[™] CMOS Sensor

Sensor				
Sensor Type	Monochrome			
Resolution	5 MP			
Live Image Resolution, No Binning	5 MP			
Live Image Frame Rate	36 f/s			
Live Image Resolution, 2 x 2 Binning	1.25 MP			
Live Image Frame Rate, 2 x 2 Binning	70 f/s			
Sensor	Sony Pregius™ CMOS, Back Illuminated Enhanced Near-IR, 6 th Generation			
Sensor P/N	Sony IMX250			
Active Pixels	2448 x 2048			
Sensor Size	8.47mm x 7.07mm			
Pixel Size	3.45µm			
Aspect Ratio	6:5			
Sensor Diagonal	11.0mm			
Optical Format	2/3"			
Electronic Shutter	Global			
Peak QE	66% at 520nm			
QE in NIR	26% at 800nm			
Camera				
Lens Mount	C-mount			
Computer Interface	USB 3.0			
Full Well	11,000 e			
Read Noise	2.1e			
Dynamic Range	74 dB			
A-D Bit Depth	12-bit			
Analog Gain	0.85, 1, 2, 4, 8			
Analog + Digital Gain	16, 32, 64, 128, 213			
Binning	2 x 2			
Cooling Technology	Thermoelectric cooler, fan-cooled heatsink, software fan/cooler disable			
Cooling	-25°C (-45°C differential from 20°C ambient)			
Temperature Readout	In SPOT software			
Dark Current	0.18 e/pixel/sec			
Exposure Time Range	28µs to 60 min			
Exposure while Downloading	Automatic			
ROI (Region of Interest)	User selectable (50 x 50 min)			

Live/Captured Image Frame Rates (frames/second)*						
Binning level	Region of Interest					
	2448 x 2048	1224 x 1024	612 x 512	306 x 256	153 x 128	
No Binning	36	69	132	242	321	
2x2 Binning	69	132	242	321	321	

*1ms manual exposure time on a computer with an Intel i5-4590s CPU (4-core processor running at 3.00 GHz).

Trigger Ready Out		
Trigger In		
Synch Out 0		
Synch Out 1		
Synch Out 2		
Synch Out 3		
(optically isolated, requires external power supply)		
Catalog: TRIG6DSUB		
External Brick Power Supply Input: AC 100-240V, 50Hz-60Hz		
Output: +5V		
20.5W		
Yes		
Catalog: QUARTER20AP		
USB 3.0 micro-B (locking)		
USB 3.0, 2m (6.5 ft), Type-A to Micro-B (locking) cable		
CE, FCC Class A, EN60950, RoHS		
Windows / Mac		
Windows / Mac / Linux		
0.55X C-Mount Adapter Catalog: DE55 Series		
RT39M5		
2.35 lbs (1.65 kg)		

Quantum Efficiency



Dimensions



.55x Adapter Magnification FOV Sony IMX250





6540 Burroughs Avenue Sterling Heights • MI • 48314-2133 phone: 586.731.6000 • fax: 586.731.6469 website: www.spotimaging.com SPOT RTCMOS Camera Data Sheet PN3173 Rev 4